




Intraoral Welder.

THE LITTLE, GREAT ALLY
OF THE IMPLANTOLOGIST.



The essential tool for implantology: it fixes and stabilizes abutments permanently, allowing prostheses to be fitted to newly placed implants in a simple, safe, and intuitive manner.


vicky/grem



Intraoral Welder.



A VERSATILE AND PRECISE TOOL.


VickyGrem's intraoral welding machine is a technologically advanced device that combines reliability and extraordinary convenience.

With the intraoral welder, it becomes fast and safe to permanently stabilize abutments connected through a titanium wire: the preset programs respond precisely to international protocols. They are customizable based on the specifications related to the individual case and the needs of the implantologist.

Easy, intuitive, and economical to use, it allows newly placed implants to be solidarized and temporary prostheses to be fitted, reconstructing total arches or partials, even in cases where there is little bone, during the same surgery.

This operation reduces patient discomfort and does not compromise osseointegration (the most delicate phase with immediate loading).

VickyGrem's intraoral welding machine is the ideal tool for safely performing immediate loading in implantology, both for one-piece implants (one-piece: welding is done directly on the abutment of the implant) as well as for implants with prosthetic components (two-pieces: welding is done on the abutment that will later be replaced with the final straight or angled abutment).



VickyGrem Intraoral Welder is an efficient instrument of high technological precision and compact size because it does not take up much space and is always ready for use. It is designed to be used intuitively, thanks partly to the touch display that guides the clinician during use.

Features:

- maximum safety: automatic detection of opening or disconnection of the clamps, control of the sequence of welding, galvanic isolation from the power supply.
- multilingual touch display: preset in Italian, English and Spanish.
- visual and audible indicators of welding steps and operations to minimize fabric heating.
- maximum titanium wire size: length 3 mm.
- reduced waiting time between one weld and the next.
- extensive customization of settings and power adjustments.

The package contains:




VG Intraoral Welder.



Ergonomic pliers autoclavable.



Foot pedal to operate without touching the display.



Connection cables, strong and easy to connect.



2-meter-long network cable.

Plus: User's Manual.

Touch display: easy and intuitive.



VickyGrem's endoral welding machine makes the implantologist's work more manageable with a touchscreen display that allows this device to be used simply and immediately. Options are presented clearly and neatly, and with a simple touch of the finger, it is possible to navigate through the welder's programs, choosing preset ones or quickly customizing parameters.

Many benefits for the implantologist and the patient.

- It increases the predictability of immediate loading.
- It lowers the risk of implant failure due to reduced implant mobility.
- It minimizes inconvenience to the patient.
- It significantly reduces time and cost.
- It allows implants to be placed and prostheses to be fitted in the same session.
- It significantly improves success rates of implant surgeries with reduced primary stability.
- Operating is simple and intuitive, thanks to the pre-set programs.
- No sparks are generated during use.
- Patient safety is total: the clamps automatically disconnect from the power supply during welding.
- The clamps can be safely autoclaved.
- The small size allows the welder to be stored away when not needed and easily transported from room to room.



With an additional advantage: VickyGrem service.

VG Intraoral Welder is not the only welder on the market, but it is the only one offered by VickyGrem, a company that does not just market products, tools, or devices. We provide solutions, support each customer based on their specific needs, complete each piece of information with advice and suggestions, and find the best way to be considered partners before suppliers.

Contact us for any questions
or information.

VickyGrem S.r.l.

www.vickygrem.com

amministrazione@vickygrem.com

tel. +39 031 879344

Via G. Mazzini, 2 -23845 - Costa Masnaga LC, Italia.

C.F. e P.IVA 03621210131 -Reg. REA LC-321965

PEC: vickygrem@legalmail.it

Sdi Recipient Code: M5UXCR1

The logo for VickyGrem features the letters 'vg' in a stylized, handwritten font above the word 'vickygrem' in a clean, lowercase sans-serif font.